



Using Illumina Genotyping Data with dChip

This document describes Illumina's data compatibility with dChip. Data sharing between BeadStudio and dChip occurs via a DLL plug-in. This plug-in allows you to use BeadStudio's Genotyping Module to create a dChip-compatible report.

For information about the dChip algorithm, visit <http://biosun1.harvard.edu/complab/dchip/>.

NOTE: Because dChip is not an Illumina product, Illumina Technical Support cannot provide support for the dChip algorithm.

Illumina Data Compatibility with dChip

Perform the following steps to import data from a dChip-compatible BeadStudio report into dChip.

1. Open dChip.
2. Get External Data.
 - a. Select **Analysis | Get External Data**.
 - b. Enter a group name.
 - c. Click **Data File** and choose the intensity and call file (e.g., projectname_dChip_raw_intensity_and_call_data.txt).
 - d. Check all three checkboxes.
 - e. Select the **Other Information** tab.
 - f. Select **External SNP data from Array Type** from the dropdown menu.
 - g. [OPTIONAL] Adjust the ploidy column in the sample info file. Use ploidy = 2 for normal samples, blank for tumor samples.
 - h. In the information files area, click in the **Sample** field and choose the sample information file (e.g. projectname_dChip_sample_info.txt).
3. Observe LOH and copy number data in dChip.
 - a. Select **Analysis | Chromosome**.
 - b. Click the **Genome Information File** field and choose the genome information file (e.g. projectname_dChip_genome_info.txt).
 - c. Select **Copy number & LOH** from the dropdown menu.
 - d. Click **Options** and adjust the settings as desired.
 - e. Click **OK** and observe LOH and copy number data in dChip.
4. [OPTIONAL] Export copy number data from dChip.
 - a. Select **Chromosome | Export SNP Data**.
 - b. Specify the output file and click **OK**.