



A. Amplicon-Based Targeted Sequencing



B. Hybridization-Based Enrichment Sequencing



C. Restriction Enzyme Reduced Representation Sequence-Based Genotyping



Figure 1: NGS Methods for Discovering and Genotyping SNPs. Amplicon-based targeted resequencing methods (panel A) adapted from Mamanova et al, 2010¹⁶ and Liu et al, 2012.¹⁸ Hybridization-based enrichment sequencing methods (panel B) adapted from Mamanova et al, 2010¹⁶ and Cronn et al, 2012.¹⁰ Restriction enzyme reduced representation sequence-based genotyping methods (panel C) adapted from Andolfatto et al, 2011.¹⁹

can be aligned, compared, and screened for SNP variants (Figure 1C).^{5,8,9,20} NGS-compatible fragment libraries enable massively parallel and multiplexed sample sequencing, facilitating the rapid discovery and genotyping of tens to hundreds of thousands of SNPs across large populations.

RE-GBS protocols, initially established for crops like maize and wheat, have advantages in cost per sample and application in species where there is no *a priori* knowledge of the genome. The application of RE-GBS is especially powerful in mapping populations, or closely related groups of samples, such as candidates for

AATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTACCCTAACGGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTAA
 AAGCTACCGTAACGAAAGTATCAATTAAGATTACTTGTATCCACTGATTCAACGTACCCTAACGGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTAA
 ACGAAAAGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTACCCTAACGGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTAA
 AGTAACGATTAAAGAGCTACCGTCAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTACCCTAACGGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTAA
 AATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTACCCTAACGGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTAA
 ATCTGATCCACTGATTCAACGTTAAGATTACTTGTATCCACTGATTCAACGTACCCTAACGGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTAA
 TATCAATTGAGACTAAATATTAACGTACCATTAAAGATTACTTGTATCCACTGATTCAACGTACCCTAACGGAACGATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCAACGTAA

genomic selection. If populations are more divergent than expected or target novel species, RE-GBS protocols can require optimization (beyond published protocols) to customize coverage and minimize missing data. For example, high divergence across targeted samples can result in missing data, complicating downstream analysis, whereas low divergence can result in a lower number of detected SNPs.

The advantages of RE-GBS are many, making the protocol development for species-specific applications rewarding.⁸ Reduced ascertainment bias over array-based methods, the ability to discover and characterize polymorphisms simultaneously, and the generation of valuable genetic information for a low (< \$20 USD) cost per sample (excluding bioinformatics) make this a method of choice for those moving from array methods to genotyping by sequencing. RE-GBS data analysis methods are supported with open-source analysis tools (eg, TASSEL) that can be tailored for crops of interest using a command-line interface. Table 2 (reproduced from Nielsen et al, 2011²¹) shows a list of available non-commercial NGS genotyping calling software. Nielsen et al also present a workflow for converting NGS data into SNP calls (Figure 2).

RAD-Seq protocol enhancements have been primarily focused on increasing the level of multiplexing to reduce cost and eliminate expensive steps in the protocol workflow, such as random shearing and the subsequent need for end repair. Examples of methods that eliminate random shearing include MSG,¹⁹ CRoPS,²² and ddRADSeq.²⁰ The ddRADSeq method has been used to refine size selection, recovering a “tunable number of regions” distributed randomly throughout the genome at a reported library preparation cost of \$5 USD per sample and input amounts as low as 100 ng of starting DNA.²⁰ This approach also implements a two-index combinatorial multiplex system ($n \times m$ individuals using $n+m$ indexes), a sequence filter analysis toolkit, and a sample tracking data management tool available through a Google Docs interface. High-throughput data management and sample tracking are critical for implementing any sample screening method in breeding and germplasm tracking.¹¹

Table 3 summarizes published sequence-based genotyping methods, including PCR-based, hybridization-based, and restriction enzyme approaches.




Figure 2: Converting NGS Data Into Genotype Calls.

Reproduced from Nielsen et al, 2011.²¹ First, pre-processing steps transform NGS data into aligned reads with quality scores that indicate confidence. Next, SNP or genotype calls are made using a multi-sample or single-sample calling procedure, depending on the number of samples and depth of coverage. Finally, post-processing steps filter the called SNPs.

AAAGAAATGATAACAGTAACACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACCGTACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAAACGACGAAAAGAATGATAACAGTAAACACACTTCTGTAAAC
AATCAACGTACCGTAACGAAACGATCAATTAAAGATTACTTGATCCACTGATTCACCGTACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAAACGACGAAAAGAATGATAACAGTAAACACACTTCTGTAAAC
AACGACGAAAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACCGTACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAAACGACGAAAAGAATGATAAC
TTAAAGTACCATTAAAGAGCTACCGTGCAAACGATTAACCAACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACCGTACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAAACGACGAAAAGAATGATAAC
AAAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACCGTACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAAACGACGAAAAGAATGATAAC
AAGATTACTTGATCCACTGATTCACCGTACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAAACGACGAAAAGAATGATAAC
AACGTTCAATTGAGACTAAATATAACGTACCATTAAAGATTACTTGATCCACTGATTCACCGTACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGCTACCGTGCAAACGACGAAAAGAATGATAAC

Table 3: Published Sequence-Based Genotyping Methods.

Method	Type of Method	Description
Amplicon sequencing	PCR based	Often used in metagenomics applications where 16S fragments are targeted. Labor intensive to amplify and tag multiple targets to optimize sequencing coverage. Difficult to scale currently to leverage sequencing output to drive down price per sample.
LR-PCR	PCR based	Long-range PCR (< 35 kbp, typically 3–10 kbp) can be used to target regions that then require shearing before library preparation. Challenges include equimolar pooling of sample/fragments. Tends to be a drop in coverage at the ends that can be resolved by increasing amplicon overlap to a minimum of 100 bp. ^{10,16}
Molecular inversion probes	Whole genome	Molecular inversion probes, single-stranded oligonucleotides with common linker flanked by target-specific sequences, anneal to target sequence and become circularized by a ligase. PCR amplification and products are sequenced directly. Suited for few targets and high sample numbers (> 100 samples). ^{25–27}
WGS/Genome skimming	Whole genome	Whole-genome sequencing includes DNA shearing and repair before adapter ligation. Low depth or genome skimming of whole genomes is performed for organelle (plastome, mitochondrial, or rDNA), phylogenetic/systematics, or comparative analysis. Can provide partial sequences of low-copy nuclear loci for designing PCR primers or probes for subsequent hybridization-based genome reduction approaches. ¹³
OS-Seq	Hybridization based	Oligonucleotide-selective sequencing is a targeted genome resequencing in which the lawn of oligonucleotide primers of an Illumina flow cell is modified to function as both a capture and sequence substrate. ²⁸
Array hybridization capture (with or without C0t1)	Hybridization based	Fragment library hybridized to immobilized probe. Non-specific hybrids are removed and targeted DNA is eluted and sequenced. Can be less labor intensive than PCR amplification. Can be followed by a target-specific array that enriches for target in a reduced-complexity sample. ^{15,16}
In-solution hybridization capture (with or without C0t1)	Hybridization based	Specific probes designed to target regions of interest from sequencing library. An excess of probes over template can result in a higher hybridization than with array-based methods. Can be more amenable to scalable throughput. ¹⁷
CRoPS	Restriction digest	Complexity reduction using AFLP with next-generation sequencing. Enables SNP discovery using tagged libraries of 2 or more genetically diverse samples. Uses a methylation-sensitive restriction enzyme sequenced at 5–10× redundancy. Use of homozygous lines is encouraged to enable selection of SNPs located in low- or single-copy genome sequences. ²²
RAD-Seq	Restriction digest	Genomic DNA digested with a restriction enzyme and a barcoded adapter is ligated to compatible sticky ends. DNA samples, each with a different barcode, are pooled, randomly sheared, and size selected (300–700 bp), and a second adapter is ligated after polishing and filling ends. A Y-adapter ensures that only RAD tags are amplified in the PCR step. ⁵
Cornell GBS	Restriction digest	Employs unmodified adapters (ie, without the 5' phosphate group and fork), removes fragment size selection. By using a single well for genomic DNA digestion and adapter ligation, it has reduced a number of enzymatic and purification steps. Methylation-sensitive enzymes are used to avoid repetitive regions of plant genomes. ⁹
Modified Cornell GBS	Restriction digest	Modifies the original Cornell GBS method by use of 2 complementary enzymes (a “rare” cutter and a “common” cutter) and a Y adapter where Adapter 1 and Adapter 2 are on opposite ends of each fragment. ⁸
ddRADSeq	Restriction digest	Relies on the concept of RAD-Seq, but eliminates the random shearing. Explicitly uses size selection to recover a tunable number of regions distributed randomly through the genome. Provides an index, computational analysis tool kit, and lightweight data management tools to facilitate multiplexing of many hundreds of individuals. Major cost reductions are attributed to removal of random shearing and subsequent end repair requirements. ²⁰
GR-RSC	Restriction digest	Genome reduction based on restriction site conservation. Includes a double digest of DNA with rare and frequent restriction enzymes, labeling a recognition rare cutter site with 5' biotin using paramagnetic bead separation, adding barcode sequences using PCR, equimolar pooling of samples, and size selection using gel isolation. ^{29–30}
MSG	Restriction digest	Multiplex NGS protocol, includes a fragment size-selection step developed to identify recombinant breakpoint of many samples simultaneously at resolution sufficient for most mapping purposes. Incorporates aspects of WGS and RAD-Seq. Uses a more frequent cutter than RAD-Seq and allows ligation of adapters to many small genomic fragments in a single step. Fragment orientation is random regarding the direction of sequencing. No shearing or repair of DNA before adapter ligation. ¹⁹
DARTSeq	Restriction digest	Based on genome complexity reduction using restriction enzymes followed by sequencing. ³¹

Summary

Genotyping arrays forged the foundation of the genomics movement in agriculture, identifying SNPs associated with desired phenotypic traits that researchers have used to improve livestock breeding and crop yields. The rapid evolution of sequencing technologies is driving the development of lower-cost sequencing-based genotyping methods that will enable agrigenomics researchers to study livestock, crops, and biological systems at a level never before possible. Providing a genome-wide view, NGG methods offer the specificity, reproducibility, and efficiency needed to accelerate agricultural research, advance the development of high-value trait screening methods, and enable the swift deployment of these applications in the real world.

References

1. Batley, J. and D. Edwards (2007). SNP applications in plants. In Oraguzie NC, Rikkerink EHA, Gardiner SE, and Silva HN (Eds.), *Association Mapping in Plants*, 95–102. New York, New York, USA: Springer.
2. Wang DG, Fan JB, Siao CJ, et al. Large-scale identification, mapping, and genotyping of single-nucleotide polymorphisms in the human genome. *Science*. 1998;280:1077-1082.
3. Van Tassel CP, Smith TPL, Matukumalli LK, et al. SNP discovery and allele frequency estimation by deep sequencing of reduced representation libraries. *Nat Methods*. 2008;5:247-252.
4. Boichard D, Chung H, Dassonneville R, et al. Design of a low-density SNP array optimized for imputation. *PLoS One*. 2012;7:e34130.
5. Baird NA, Etter PD, Atwood TS, et al. Rapid SNP discovery and genetic mapping using sequenced RAD markers. *PLoS One*. 2008;3:e3376.
6. Kirst M, Resende M, Munoz P and Neves L. Capturing and genotyping the genome-wide genetic diversity of trees for association mapping and genomic selection. *BMC Proceedings*. 2011;5:17.
7. Metzger ML. Sequencing technologies—the next generation. *Nat Rev Genet*. 2010;11:31-46.
8. Poland JA, Brown PJ, Sorrells ME, Jannink J. Development of high-density genetic maps for barley and wheat using a novel two-enzyme genotyping-by-sequencing approach. *PLoS One*. 2012;7:e32253.
9. Elshire RJ, Glaubitz JC, Sun Q, et al. A robust, simple genotyping-by-sequencing (GBS) approach for high diversity species. *PLoS One*. 2011;6:e19379.
10. Cronn R, Knaus BJ, Liston A, et al. Targeted enrichment strategies for next generation plant biology. *Am J Bot*. 2012;99:291-311.
11. McCouch SR, McNally KL, Wang W, Sackville Hamilton R. Genomics of gene banks: A case study in rice. *Am J Bot*. 2012;99:407-423.
12. Lorenc MT, Hayashi S, Stiller J, et al. Discovery of single nucleotide polymorphisms in complex genomes using SGSautoSNP. *Biology*. 2012;1:370-382.
13. Huang X, Feng Q, Quian Q, et al. High-throughput genotyping by whole-genome resequencing. *Genome Res*. 2009;19:1068-1076.
14. Grover CE, Salmon A, Wendel JF. Targeted sequence capture as a powerful tool for evolutionary analysis. *Am J Bot*. 2012;99:312-319.
15. Fu Y, Springer NM, Gerhardt DJ, et al. Repeat subtraction-mediated sequence capture from a complex genome. *Plant J*. 2010;62:898-909.
16. Mamanova L, Coffey AJ, Scott CE, et al. Target-enrichment strategies for next-generation sequencing. *Nat Methods*. 2010;7:111-118.
17. Saintenac C, Jiang D, Akhunov ED. Targeted analysis of nucleotide and copy number variation by exon capture in allotetraploid wheat genome. *Genome Biol*. 2011;12:R88.
18. Liu S, Yeh CT, Tang HM, Nettleton D, Schnable PS. Gene mapping via bulked segregant RNA-Seq (BSR-Seq). *PLoS One*. 2012;7:e36406.
19. Andolfatto P, Davison D, Erezylmaz D, et al. Multiplexed shotgun genotyping for rapid and efficient genetic mapping. *Genome Res*. 2011;21: 610-617.
20. Peterson BK, Weber JN, Kay EH, Fisher HS, Hoekstra HE. Double digest RADSeq: an inexpensive method for *de novo* SNP discovery and genotyping in model and non-model species. *PLoS One*. 2012;7:e37135.
21. Nielsen R, Paul JS, Albrechtsen A, Song YS. Genotype and SNP calling from next-generation sequencing data. *Nat Rev Genet*. 2011;12:443-451.
22. Van Orsouw NJ, Hogers RC, Janssen A, et al. Complexity reduction of polymorphic sequences (CRoPS): a novel approach for large-scale polymorphism discovery in complex genomes. *PLoS One*. 2007;2:e1172.
23. Li Y, Sidore C, Kang HM, Boehnke M, Abecasis GR. Low-coverage sequencing: implications for design of complex trait association studies. *Genome Res*. 2011;21:940-951.
24. 1000 Bull Genomes Project (www.1000bullgenomes.com) Accessed 18 December 2014.
25. Hardenbol P, Baner J, Jain M, et al. Multiplexed genotyping with sequence-tagged molecular inversion probes. *Nat Biotechnol*. 2003;21:673-678.
26. Hardenbol P, Yu F, Belmont J, et al. Highly multiplexed molecular inversion probe genotyping: over 10,000 targeted SNPs genotyped in a single tube assay. *Genome Res*. 2005;15:269-275.
27. Dahl F, Gullberg M, Stenberg J, Landegren U, Nilsson M. Multiplex amplification enabled by selective circularization of large sets of genomic DNA fragments. *Nucleic Acids Res*. 2005;33:e71.
28. Mylykangas S, Buenrostro JD, Natsoulis G, Bell JM, Ji HP. Efficient targeted resequencing of human germline and cancer genomes by oligonucleotide-selective sequencing. *Nat Biotechnol*. 2011;29:1024-1029.
29. Maughan PJ, Yourstone SM, Jellen EN, Udall JA. SNP discovery via genomic reduction, barcoding and 454-pyrosequencing in *Amaranth*. *Plant Genome*. 2009;2:260-270.

30. Maughan PJ, Yourstone SM, Byers RL, Smith SM, Udall JA. Single-nucleotide polymorphism genotyping in mapping populations via genomic reduction and next-generation sequencing: proof of concept. *Plant Genome*. 2010;3:166-178.
31. Sansaloni C, Petroli C, Jaccoud D, et al. Diversity arrays technology (DArT) and next-generation sequencing combined: genome-wide, high throughput, highly informative genotyping for molecular breeding of Eucalyptus. *BMC Proceedings*. 2011;5:P54.
32. Gore MA, Wright MH, Ersoz ES, et al. Large-scale discovery of gene-enriched SNPs. *Plant Genome*. 2009;2:121-133.
33. Deschamps S, la Rota M, Ratashak JP, et al. Rapid genome-wide single nucleotide polymorphism discovery in soybean and rice via deep resequencing of reduced representation libraries with the Illumina Genome Analyzer. *Plant Genome*. 2010;3:53-68.
34. Straub SCK, Parks M, Weitemier K, Fishbein M, Cronn RC, Liston A. Navigating the tip of the genomic iceberg: next-generation sequencing for plant systematics. *Am J Bot*. 2012; 99:349-364.

AAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA
AATCAACGTACCGTAAACGAACGTATCAATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA
AACGACGAAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA
TTAAGGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA
AAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA
AAGATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA
AACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA

