

## Quantify and Dilute DNA

- 1 Quantify DNA using a fluorometric method.
- 2 Dilute DNA to 25 ng/μl in RS1.
- 3 Requantify the diluted DNA.
- 4 In a microcentrifuge tube, dilute the 25 ng/μl DNA in RS1 to a total volume of 4 μl.
- 5 Add 1 μl SS1 to the 4 μl of 50 ng DNA.

## Hybridize Oligo Pool

- 1 Dilute 2.5 μl CAT with 2.5 μl RS1 per sample well.
- 2 Pulse vortex to mix, and then centrifuge briefly.
- 3 Add 5 μl RS1 to one well.
- 4 Add 5 μl diluted DNA to the remaining wells.
- 5 Add 5 μl diluted CAT to all wells.
- 6 Add 15 μl OHS2 to each well. Pipette to mix.
- 7 If bubbles form, centrifuge the plate at 100 × g for 20 seconds.
- 8 Place on the thermal cycler and run the HYB program.
- 9 For 96 samples, combine ELE and ELB as follows.
  - a Transfer 137 μl ELE to the contents of the ELB tube.
  - b Flick and invert to mix.
- 10 Place the ELB/ELE mixture on ice for later use.

## Remove Unbound Oligos

- 1 Add 25 μl SPB. Pipette slowly to mix.
- 2 Incubate at room temperature for 5 minutes.
- 3 Place on the magnetic stand until liquid is clear.
- 4 Remove and discard all supernatant.
- 5 Wash three times with 80 μl SW1.
- 6 Use a 20 μl pipette to remove residual SW1.
- 7 Add 80 μl of 60% EtOH.
- 8 Incubate at room temperature for 30 seconds.
- 9 Remove and discard all supernatant.
- 10 Use a 20 μl pipette to remove residual EtOH.
- 11 Air-dry for up to 5 minutes.

## Extend and Ligate Bound Oligos

- 1 Remove the plate from the magnetic stand.
- 2 Add 22 µl ELB/ELE mixture to each well.
- 3 Pipette to mix.
- 4 If bubbles form, centrifuge at 100 × g for 20 seconds.
- 5 Place on the thermal cycler and run the EXT\_LIG program.
- 6 Combine EDP and EMM as follows.
  - ▶ [1 sample] 1.1 µl EDP and 21 µl EMM
  - ▶ [96 samples] 106 µl EDP and 2006 µl EMM
- 7 Pipette the EDP/EMM mixture to mix, and then centrifuge briefly.
- 8 Place the EDP/EMM mixture on ice for later use.

## Amplify Libraries

- 1 Arrange the Index 1 adapters in columns 1–12.
- 2 Arrange the Index 2 adapters in rows A–H.
- 3 Place the HYP plate on a TruSeq Index Plate Fixture.
- 4 Add 4 µl of each Index 1 adapter down each column.
- 5 Add 4 µl of each Index 2 adapter across each row.
- 6 Place the plate on ice or iceless cooler.
- 7 Add 20 µl EDP/EMM mixture. Pipette to mix.
- 8 Centrifuge at 280 × g for 1 minute.
- 9 Place the plate on ice or iceless cooler.
- 10 Immediately transfer to the post-PCR area.
- 11 Place on the thermal cycler and run the PCR program for the appropriate number of cycles.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

## Clean Up Libraries

- 1 Centrifuge the HYP plate at 280 × g for 1 minute.
- 2 Add 36 µl SPB to the CLP plate.
- 3 Place the HYP plate on a magnetic stand until liquid is clear.
- 4 Transfer 45 µl supernatant from the HYP plate to the CLP plate.
- 5 Shake the plate at 1800 rpm for 2 minutes.
- 6 Incubate at room temperature for 5 minutes.
- 7 Centrifuge at 280 × g for 1 minute.
- 8 Place on a magnetic stand until liquid is clear.
- 9 Remove and discard all supernatant.
- 10 Wash two times with 200 µl 80% EtOH.
- 11 Using a 20 µl pipette, remove residual EtOH.
- 12 Remove from the magnetic stand and air-dry for up to 5 minutes.
- 13 Add 25 µl RSB.
- 14 Shake the plate at 1800 rpm for 2 minutes.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge at 280 × g for 1 minute.
- 17 Place on a magnetic stand until liquid is clear.
- 18 Transfer 20 µl purified library from the CLP plate to the LNP plate.
- 19 From the liquid in the CLP plate, run an aliquot of the samples and control to confirm the PCR product sizes.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 6 months.

## Normalize Libraries

- 1 Add 44 µl LNA1 per library to a new 15 ml conical tube.
- 2 Use a P1000 pipette to resuspend LNB1.
- 3 Transfer 8 µl LNB1 per library to the tube of LNA1. Invert to mix.
- 4 Add 45 µl LNA1/LNB1 to the LNP plate.
- 5 Shake at 1800 rpm for 30 minutes.
- 6 Place on a magnetic stand until liquid is clear.
- 7 Remove and discard all supernatant.
- 8 Remove from the magnetic stand.
- 9 Wash two times with 45 µl LNW1. Shake at 1800 rpm for 5 minutes per wash.
- 10 Use a 20 µl pipette to remove residual LNW1.
- 11 Remove from the magnetic stand.
- 12 Add 30 µl fresh 0.1 N NaOH.
- 13 Shake at 1800 rpm for 5 minutes.
- 14 Place on a magnetic stand until liquid is clear.
- 15 Add 30 µl LNS2 to the SGP plate.
- 16 Transfer 30 µl supernatant from the LNP plate to the SGP plate.
- 17 Centrifuge at 1000 × g for 1 minute.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 30 days.

## Pool Libraries

- 1 Centrifuge at 1000 × g for 1 minute.
- 2 Transfer 5 µl of each library to an 8-tube strip.
- 3 Seal the plate and store at -25°C to -15°C.
- 4 Transfer the contents of the 8-tube strip to the PAL tube. Pipette to mix.
- 5 Denature and dilute the library pool to the appropriate loading concentration.

### SAFE STOPPING POINT

If you are stopping, cap the tubes and store at -25°C to -15°C for up to 7 days.

## Acronyms

| Acronym | Definition                                   |
|---------|--|
| CLP     | Cleanup Plate                                |
| EDP     | Enhanced DNA Polymerase                      |
| ELB     | Extension-Ligation Buffer                    |
| ELE     | Extension-Ligation Enzyme                    |
| EMM     | Enhanced Master Mix                          |
| CAT     | TruSeq Genotype N <sub>e</sub> Oligos        |
| HT1     | Hybridization Buffer                         |
| HYP     | Hybridization Plate                          |
| LNA1    | Library Normalization Additives 1            |
| LNB1    | Library Normalization Beads 1                |
| LNP     | Library Normalization Plate                  |
| LNS2    | Library Normalization Storage Buffer 2       |
| LNW1    | Library Normalization Wash 1                 |
| OHS2    | Oligo Hybridization for Sequencing Reagent 2 |
| PAL     | Pooled Amplicon Library                      |
| RS1     | Resuspension Solution 1                      |
| RSB     | Resuspension Buffer                          |
| SGP     | Storage Plate                                |
| SNP     | Single Nucleotide Polymorphism               |
| SPB     | Sample Purification Beads                    |
| SS1     | Sample Stabilization Solution 1              |
| SW1     | Stringent Wash 1                             |