

DRAGEN TSO500 ctDNA Analysis Software on ICA

Release Notes

V2.5.0

For TruSight Oncology 500 ctDNA Assay

March 11, 2024 March 6, 2024

Introduction

These Release Notes detail the key features and known limitations of software components for the DRAGEN TSO500 ctDNA v2.5.0 Analysis Software on ICA. Below is a summary of the changes included in DRAGEN TSO 500 ctDNA v2.5.0 Software on ICA. For full details, please consult the DRAGEN TSO 500 ctDNA v2.5.0 on ICA Software User Guide available on the support website.

This software is intended for use with the TruSight Oncology 500 ctDNA and TruSight Oncology 500 ctDNA v2 assays.

- Software Version: 2.5.0
- DRAGEN software version 3.10.15

NEW FEATURES:

ICA functionality

- Enabled per-sample-per-node processes which reduce the turnaround time that was observed when multiple long running samples were grouped in the same node.
- Enabled easy access to NextFlow process logs to facilitate troubleshooting.

BaseSpace Sequence Hub (BaseSpace) functionality

- No new functionality

Bioinformatics pipeline

- All SNVs, insertions and deletions that are part of MNVs are now reported both individually and as merged variants in the final VCF and the Combined Output Variant file.
- DRAGEN sex prediction algorithm was added to the TSO500 pipeline. The sex is predicted based on the read count information in the sex chromosomes and the autosomal chromosomes.
- MSI JSON file containing MSI results was added to the Results folder.
- SV evidence BAMs were added to Logs_Intermediates folder.

FIXED ISSUES:

- Fixed an out of memory issue that was impacting samples that were sequenced at higher (than recommended) coverage on AWS cloud instances.
- Illumina Annotation Engine 3.2.6 (aka Nirvana) includes the following bug fix:
 - A fixed RefSeq version (105.20220307) was incorporated that fixed canonical transcript assignments for some prominent genes and variants. For example, for variant BRAF NP_004324.2. V600E, the canonical transcript and HGVS notations are now fixed in the CombinedVariantOutput file (was presented as BRAF NP_001361187.1. V640E in the previous version).
- Fixed an error of CombinedVariantOutput file displaying more than one annotation for some variants.
- Fixed errors related to the sample sheet validator in BaseSpace Run Planning tool.
- Fixed an error of analysis failure when the input folder name has only numeric characters.
- Fixed an error of DRAGEN logs not being saved in the output directory.
- Fixed an error in the MetricsOutput file where NTC samples (No-Template Control samples with 0 reads) were indicated as "Pass" instead of "Fail".
- Fixed errors with MNV calling when phased variant merging distance is greater than 10bp. Ensure haplotype from sufficiently long k-mers are used for merging phased variants.
- Fixed an instance where DRAGEN small variant calling did not correctly handle overlapping mates, mistakenly detected strand bias and subsequently incorrectly filtered the variant. Strand bias was measured because the strongest support for the variant consistently came from one strand. This should however not be measured as strand bias since both mates were in agreement.
- Fixed an issue causing runs with samples with extreme copy number gains (e.g. fold change > 50, corresponding to ~250 copies when tumor fraction is ~40%) in a particular region or contrived samples to take significantly longer than 20 hours.
- Updated SV caller for more efficient SV processing.

KNOWN ISSUES:

- The sample sheet should not have blank rows between samples in the [Data] section, this may cause a run failure.
- Performance not verified using reads other than 2 x 151, paired end, dual index.
- The software does not notify the user when InterOp files for RunQC are missing or corrupted.
- Some contrived samples such as SeraCare Complete Mutation Mix, which have multiple structural variants (SVs) and high library conversion efficiencies, could generate a high number of chimeric reads and high number of candidate SVs. Occasionally, the SV caller may filter some of the reads and lead to occasionally missing fusions. In such cases downsampling the FASTQs can help recover those fusion calls. Contact your local support team for additional details and a workaround.

- ICA does not provide sufficient debugging information, when the user auto-launches the TSO 500 ctDNA v2.5 pipeline and the input run folder contains additional, unexpected data that reduces available space for analysis.
- Analysis fails when starting from V1 sample sheets due to missing adapter sequences in V1 sample sheet template. Users are recommended to start with V2 sample sheet template or add adapter sequences manually.
- Pipeline does not exit early and continues to the next DragenCaller step due to TSO500 ctDNA FASTQ validation failure if Fastq_list.csv is missing.
- When storage is not selected as "Large", FastqGeneration step fails when using S4 flow cells and the run folder size is up to 1.2 TB, and FASTQ files up to 4.2 TB. This is due to FASTQ files being duplicated in both the Nextflow works folder as well as the Logs_Intermediates/FastqGeneration folder causing the disk space to run out before the FastqGeneration step could be completed.
- ctDNA pipeline fails for an NTC sample (No-Template Control samples with 0 reads) due to absence of Evidence BAM File.
- In the V2 CNV cutoff bed file, gene "MYCL" should be listed instead of "MCYL1".

PRODUCT LIMITATIONS:

- The sample sheet must be configured as described in [the provided templates](#), the User Guide or by using BaseSpace Run Planning tool.
- Sample sheets generated for auto-launch on ICA are not compatible and cannot be reused without changes for DRAGEN TSO500 ctDNA Analysis Software on a Local DRAGEN server, and vice versa.
- ICA run time depends on ICA instance availability, it will be affected by region and traffic
- Added validation for the storage size selection generates an error if "Small" or "Medium" values are selected ("Large" is required as a minimum) but the error message appears with a delay.
- The values in the Run Metrics section will be listed as 'NA' if the analysis was started from FASTQs or if the analysis was started from BCLs but the InterOp files are missing or corrupted.
- Germline estimation uses the latest publicly available population data and is estimated to be representative of targeted population, the impact of rare germline mutations is expected to be limited.
- The Illumina Annotation Engine (aka Nirvana) may report incorrect HGVS c. and HGVS p. notation for small variants occurring in RefSeq transcripts that exhibit transcript sequences differing from the genomic reference (i.e., RNA-edits). Currently the HGVS c. error rate is 0.00527% and the HGVS p. error rate is 0.00737%.
- The CNV caller has slightly higher noise for sample types that are not included in the baseline used for normalization (eg., cell lines). The baseline samples consist of mostly healthy donor clinical samples and SeraCare-contrived samples.
- MSAF output has had limited testing and needs to be used with caution. Updates to the small variant calling have led to an increased MSAF in samples with higher DNA input.

Release History

Version	Workflow#	Author	Description of Change
00	CN 1104006	Svetlana Bureeva	Initial Release