

# Maximizing performance on the MiSeq™ i100 Series

Library loading optimization  
steps to ensure run success

- Determine optimal library loading concentration for MiSeq i100 Series flow cells
- Improve insert size representation to maximize sequencing performance
- Support low-diversity libraries by adjusting library complexity with PhiX



## Introduction


The MiSeq i100 Series offers the simplest, fastest benchtop sequencing. Breakthrough advancements in system design, XLEAP-SBS™ chemistry, and integrated data analysis deliver enhanced usability, high data accuracy, and exceptional speed, generating results up to 4× faster than the original MiSeq System. As part of an end-to-end NGS solution, the MiSeq i100 Series provides same-day results for various applications, including transcriptomics, microbial genomics, and targeted gene sequencing studies in key areas, such as microbiology, infectious disease, oncology, and more.

When transitioning projects to the MiSeq i100 Series from another sequencing system, optimizing library loading can help maximize data yield and quality. This technical note provides recommendations to optimize results on the MiSeq i100 Series, including guidance on library loading concentration, library quality, and nucleotide diversity considerations.

## Optimal library loading

Loading concentration refers to the final concentration of a library loaded onto an instrument for sequencing. After libraries are prepared, they are diluted to the loading concentration appropriate for the library type, sequencing system, and reagent kit.

Loading libraries at a concentration that is too high or too low may lead to lower sequencing quality and yield, and, possibly, run failures in extreme conditions. Underloading may result in a low percentage of nanowell occupancy (% Occupied) and higher duplicate reads, which then requires more reads to achieve target coverage. In contrast, overloading may result in a low percentage of clusters passing filter (% PF). To determine the optimal loading concentrations on the MiSeq i100 Series, the % Occupied and % PF metrics can be plotted in Sequencing Analysis View to determine if a run was underloaded, optimally loaded, or overloaded. The approach in the following example experiment can be used to titrate the loading concentration to assess primary and secondary metrics.

 To learn more, read [Optimizing library loading for Illumina NGS systems with patterned flow cells](#)

## Determining optimal loading concentration

When finding the optimal loading concentration, it is critical to test a wide range of concentrations. Use primary metrics like % PF and % Occupied with secondary metrics like duplicates, insert size, and coverage to measure performance at various loading concentrations to determine the "usable yield" for a given application.

### Step 1: Design titration experiment

For transitioning projects from the original MiSeq System using the MiSeq Reagent Kit v3 to the MiSeq i100 Series, center titrations at ~6.5× the MiSeq Reagent Kit v3 loading concentration. Recommended centerpoint concentrations vary for different library preparation kits for use with the MiSeq i100 Series ([Table 1](#)). For all other cases, it is recommended to use 100 pM for the centerpoint concentration.

**Table 1: Recommended centerpoint concentrations for titration design with the MiSeq i100 Series<sup>a</sup>**

Library preparation kit	Centerpoint concentration
Illumina DNA Prep	80 pM
Illumina DNA Prep with Enrichment	60 pM
Illumina RNA Prep with Enrichment	80 pM
Illumina DNA PCR-Free	120 pM
TruSeq™ DNA PCR-Free	120 pM
TruSeq DNA Nano	120 pM
Illumina Viral Surveillance Panel v2	80 pM
Illumina Microbial Amplicon Prep—Influenza A/B	80 pM
Respiratory Pathogen ID/AMR Enrichment Panel	80 pM
Urinary Pathogen ID/AMR Panel	80 pM
TruSight™ RNA Pan-Cancer	80 pM
16S rRNA Amplicon <sup>b</sup>	60 pM
Pillar® oncoReveal™ Myeloid Panel	80 pM
Pillar oncoReveal Essential MPN Panel	80 pM
Pillar oncoReveal Multi-Cancer v4 with CNV Panel	80 pM
Pillar oncoReveal BRCA1 & BRCA2 plus CNV Panel	80 pM
PhiX Control v3	120 pM

- a. Double-stranded DNA libraries were quantified using the fluorometric Qubit dsDNA Quantitation High Sensitivity assay (Thermo Fisher, Catalog no. Q32851), and Bioanalyzer High Sensitivity DNA Kit (Agilent, Catalog no. 5067-4626) for average fragment size estimation. Single-stranded DNA libraries were quantified using the Qubit ssDNA Assay Kit (Thermo Fisher, Catalog no. Q10212).
- b. 16S rRNA Amplicon libraries prepared using workflow described in document [16S Metagenomic Sequencing Library Preparation](#) (Part # 15044223 Rev.B).

In this example, a library pool consisting of bacterial genome samples from *Bacillus pacificus*, *Cereibacter sphaeroides*, and *Escherichia coli* prepared using Illumina DNA Prep was tested at loading concentrations of 40 pM, 80 pM, and 120 pM.

#### Step 2: Assess nanowell occupancy and clusters PF

Plot the % PF vs % Occupied metrics from the sequencing run for each loading concentration to determine which concentrations resulted in underloading, overloading, or balanced loading. In this example: all three concentrations tested (40 pM, 80 pM, 120 pM) display optimal loading shape (a cloud of points with a positive slope) in the % PF vs % Occupied plot, demonstrating that the MiSeq i100 Series can achieve robust results within a broad library loading concentration range (Figure 1).

#### Step 3: Assess duplicates

Narrow the target concentration range by analyzing the percent of duplicates. Duplicates tend to decrease with increasing loading concentration. In this example, while all three concentrations tested have duplicates less than 15%, 80 pM and 120 pM had the lowest amount (Figure 2).

#### Step 4: Analyze insert size

Review the insert sizes. The optimal range for your library and application may vary depending on your workflow requirements. In this example, insert sizes for all three bacterial strains vary across the concentration range tested, with the greatest difference observed between 40 pM and 80 pM (Figure 2).

#### Step 5: Review other application-dependent metrics (coverage, mapping, etc)

Review additional secondary analysis metrics for optimal performance of your application. In this example, the percent-mapped metric shows robust results for all three loading concentrations tested (Figure 2). Secondary metrics were generated with the DRAGEN™ Small Whole Genome Sequencing app that is available as an on-instrument and on-cloud solution.




Figure 1: Optimal nanowell occupancy across a broad library loading concentration range—Sequencing of libraries loaded at 40 pM, 80 pM, and 120 pM displayed an optimal loading shape, demonstrating that the MiSeq i100 Series can achieve robust results across a broad range of loading concentrations.




Figure 2: Optimizing sequencing on the MiSeq i100 Series—Example titration experiment looking at duplicates, average coverage, and insert size.

## Library quality

Short inserts and contaminants introduced during library preparation, including adapter dimers, primer dimers, and partial library constructs, can negatively impact clustering on the MiSeq i100 Series. It is critical to remove these short inserts and contaminants during cleanup or size-selection steps. If necessary, short inserts and contaminants can be more effectively removed by adding an optional bead purification step to the library preparation protocol. After

library preparation is complete and before sequencing, users should verify the quality and purity of all libraries. Use an Agilent Bioanalyzer, Fragment Analyzer system, or TapeStation to check for library integrity, average insert size, and contaminants.

## Removing short inserts with additional bead purification

In this example, libraries of wastewater samples prepared with the Illumina Viral Surveillance Panel v2 Kit were treated with an additional round of bead purification using a bead to sample ratio of 0.8× following a protocol similar to that described in the [Optimal variant calling with Illumina DNA PCR-Free Prep on the NovaSeq™ X Series technical note](#), the only deviation being in step 2, where 40 µl of Illumina Purification Beads were added, rather than the specified 52.5 µl. The additional round of bead purification effectively removed most of fragments < 250 bp, with a total reduction in library yield of ~ 35% (Figure 3).

## Improving sequencing performance

Viral Surveillance Panel v2 libraries, with and without the additional round of bead purification, were sequenced on the MiSeq i100 Series and analyzed with the DRAGEN™ Microbial Enrichment Plus app. Sequencing of libraries prepared with the additional round of bead purification resulted in improved metrics, including mean read length and % post-quality reads, and increased microorganism detection, as compared to the unmodified protocol (Figure 4).

## Nucleotide diversity

Nucleotide diversity indicates the relative proportion of each base (A, C, G, or T) present in every cycle of the run. Nucleotide balance is important for color matrix correction and intensity normalization by the sequencing system. The adaptive Real-Time Analysis software onboard the MiSeq i100 Series has been carefully developed for accurate basecalling of low-diversity libraries. Optimal performance of low-diversity library sequencing can be achieved with a minimum % PhiX spike-in ( $\geq 5\%$ ) to maximize the number of high-quality reads.

In this example, low-diversity 16S amplicon libraries with 5% and 20% PhiX spike-in sequenced on the MiSeq i100 Series show robust performance comparable to the performance of high-diversity human Illumina DNA Prep libraries (Figure 5).




Figure 3: Increased insert size with additional bead purification— An additional round of bead purification (red line) effectively removed most fragments less than 250 bp, compared to the unmodified protocol (blue line).




Figure 4: Improved performance with increased insert sizes— Sequencing of libraries generated following a modified protocol (with increased insert sizes) on the MiSeq i100 Series resulted in improved performance, including increased mean read length, % post-quality reads, and microorganism detection.




Figure 5: Support for low-diversity libraries—Software onboard the MiSeq i100 Series optimizes sequencing performance for low-diversity libraries, as seen by % bases  $\geq$  Q30 and Gb output. All runs were sequenced at  $2 \times 301$  bp read length using the MiSeq i100 Series 25M Reagent Kit (600 cycles), with dashed lines representing performance specifications.

## Summary

Breakthrough advancements in sequencing chemistry and integrated data analysis on the MiSeq i100 Series deliver enhanced usability, high data accuracy, and exceptional speed. Following the best practices outlined in this technical note to assess library quality, optimize loading concentration, and pool libraries can maximize performance on the MiSeq i100 Series.

## Learn more

[MiSeq i100 and MiSeq i100 Plus Sequencing Systems](#)



1.800.809.4566 toll-free (US) | +1.858.202.4566 tel  
techsupport@illumina.com | www.illumina.com

© 2024 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners.  
For specific trademark information, see [www.illumina.com/company/legal.html](http://www.illumina.com/company/legal.html).

M-GL-03322 v1.0