

PGx Cloud Analysis Release Notes

1.1.1

July 2021

Introduction

These Release Notes detail the key changes for the Arrays Analysis Workflows (ACWF).

The PGx analysis is only intended for use with the Global Diversity Array with Enhanced PGx Content v1.0 product.

NEW FEATURES:

- Standard configurations can be used in launching analyses.

DEFECT REPAIRS:

- None

KNOWN ISSUES:

- Minor discordance may be observed with certain samples between genome build GRCh 37 and GRCh 38.